



# Rollator Guide

---

Brant Arts is an ADP Funding Authorized Vendor

Please Contact Our Experts  
With Any Questions:

**1-888-678-3144**

[mobility@brantarts.ca](mailto:mobility@brantarts.ca)  
[bamobility.ca](http://bamobility.ca)



# Category 3 Rollators:

## Trillium



Size	Weight	Seat height	Wheel size	Max user weight	Seat width	Colour selection	Handle Height
Mini	14 lbs	18"	8"	300 lbs	18"	Black, Blue, Pink	29-32"
Regular	14 lbs	21"	8"	300 lbs	18"		31.5-34.5"
Tall	14 lbs	24"	8"	350 lbs	18"		34.5-37.5"

Notes: Cable free, offset wheels for more foot space, angled backrest that supports from lower back, padded leatherette seat cushion

## Nexus 3



Size	Weight	Seat height	Wheel size	Max user weight	Seat width	Colour selection	Handle Height
Super low	16.9 lbs	19"	8"	441 lbs	17.5"	Black, Blue, Red	30.5-34.1"
Low	17.1 lbs	22"	8"	441 lbs	17.5"		30.5-36.5"
Standard	17.3 lbs	24"	8"	441 lbs	17.5"		32.5-39.8"

Notes: Seat is a bit firmer and measures taller. Brakes are a bit firmer, Offset wheels for more foot space, Post style hand grips allow for easier transfers

## Xpresso



Size	Weight	Seat height	Wheel size	Max user weight	Seat width	Colour selection	Handle Height
Mini	15 lbs	18"	8"	350 lbs	18"	Black, red, blue	30.5-33.5"
Regular	15 lbs	21"	8"	350 lbs	18"		30.5-34.5"
Tall	15 lbs	24"	8"	350 lbs	18"		33.5-40"

Notes: Cable Free, Smaller footprint, Open removable Basket, Thin padded backstrap

## Stella



Size	Weight	Seat height	Wheel size	Max user weight	Seat width	Colour selection	Handle Height
18	14 lbs	18"	8"	300 lbs	18"	Black, blue, red	19.5-33.5"
20	14 lbs	20"	8"	300 lbs	18"		30.5-34.5"
23	14 lbs	23"	8"	300 lbs	18"		33.5-37.5"

Notes: Wide supportive back strap, Post handles for easier transfers, No frame near knees for more open feel when seated, Sling seat sits a little lower and overhangs basket. Comes with a cup holder and cane holder. Locks closed when folded, Built in curb climber.

## Nitro



Size	Weight	Seat height	Wheel size	Max user weight	Seat width	Colour selection	Handle Height
Hemi	17.5 lbs	18"	8"+10"	300 lbs	18"	Black, white, red	28-31"
Regular	17.5 lbs	20.5"	8"+10"	300 lbs	18"		33.5-38.5"
Tall	17.5 lbs	23.6"	8"+10"	300 lbs	18"		36-41"

Notes: Post handles for easier transfers, Sling seat sits a little lower than measured, back strap is adjustable in height. Larger front wheels for curbs, Has a locking clip, Built in curb climber.

## Excursion



Size	Weight	Seat height	Wheel size	Max user weight	Seat width	Colour selection	Handle Height
x18	15.5 lbs	18"	8"	300 lbs	17.5"	Black, red	28-31"
x20	15.5 lbs	20"	8"	300 lbs	17.5"	Red, purple, black, red	33.5-38.5"
x23	15.5 lbs	23"	8"	300 lbs	17.5"	Blue, Red, Silver, Black	36-41"

Notes: Sling seat measures a little lower, Has cables that can catch, Ergo comfort hand grips, Built in curb climber

## Convertible Options

### Stella



**(SIGNIFICANT UPCHARGE FOR CONVERTIBLE ACCESSORIES + 6lbs)**

Size	Weight	Seat height	Wheel size	Max user weight	Seat width	Colour selection	Handle Height
18	14 lbs	18"	8"	300 lbs	18"	Black, Blue, Red	19.5-33.5"
20	14 lbs	20"	8"	300 lbs	18"		30.5-34.5"
23	14 lbs	23"	8"	300 lbs	18"		33.5-37.5"

Notes: Wide supportive back strap, Post handles for easier transfers, No frame near knees for more open feel when seated, Sling seat sits a little lower and overhangs basket. Comes with a cup holder and cane holder. Locks closed when folded, Built in curb climber

### Fusion 2 in 1



Size	Weight	Seat height	Wheel size	Max user weight	Seat width	Colour selection	Handle Height
F18	17.1 lbs	18"	8"	300 lbs	19.5"	Cranberry	29.5-35"
F20	17.1 lbs	20"	8"	300 lbs	19.5"	Blue, Black, cranberry	30-36.5"
F23	19.9 lbs	23"	8"	300 lbs	19.5 "	Black, Blue	33.5-39.5 "

Notes: Built in curb climber, Exterior cables that can catch on items, **wider seat**, ergo handles for comfort when pushing, Wide backrest pad, locking clip, basket tucked under seat a bit.

## Heavy Duty Options

### Stella



Size	Weight	Seat height	Wheel size	Max user weight	Seat width	Colour selection	Handle Height
20	14 lbs	20"	8"	400 lbs	19.5"	Black	30.5-35"

Notes: Wide supportive back strap, Post handles for easier transfers, No frame near knees for more open feel when seated, Sling seat sits a little lower and overhangs basket. Comes with a cup holder and cane holder. Locks closed when folded, Built in curb climber

### Nitro HD



Size	Weight	Seat height	Wheel size	Max user weight	Seat width	Colour selection	Handle Height
Nitro HD	23 lbs	18"	8"+10"	450 lbs	23"	Black	36-41"

Notes: Wide and tall, sling seat measures a little lower, adjustable height backrest, shocks for smoother ride and high weight capacity.

### Xpresso Wide



Size	Weight	Seat height	Wheel size	Max user weight	Seat width	Colour selection	Handle Height
Low Wide	30 lbs	19"	8"	450 lbs	19.5"	Black, Red, Silver	30.5-37.5"
Wide	31 lbs	22"	8"	450 lbs	19.5"	Black, Red, Silver	32.5-37.5"

Notes: Steel reinforced frame, Heavy, Wider frame makes maneuvering a little difficult

### Excursion XWD



Size	Weight	Seat height	Wheel size	Max user weight	Seat width	Colour selection	Handle Height
XWD	16.4 lbs	20"	8"	400 lbs	19.5"	Black, Blue	30-36.5"

Notes: Sling seat measures a little lower, Has cables that can catch, Ergo comfort hand grips, Built in curb climber.

# Category 2 Rollators:

## Trillium 6



Size	Weight	Seat height	Wheel size	Max user weight	Seat width	Colour selection	Handle Height
Mini	13 lbs	18"	6"	300 lbs	18"	Black	29-32"
Regular	13 lbs	21"	6"	300 lbs	18"		31.5-34.5"
Tall	14 lbs	24"	6"	300 lbs	18"		34.5-37.5"

Notes: Wider wheels, Thin backstrap, Basic seat with vinyl cover, Open basket

## Nexus One



Size	Weight	Seat height	Wheel size	Max user weight	Seat width	Colour selection	Handle Height
17	17 lbs	17"	6"	353 lbs	17.3"	Black	28.4-31.5"
20	17.4 lbs	20"	6"	353 lbs	17.3"		29-35"
22	17.9 lbs	22"	6"	353 lbs	17.3"		31-38"

Notes: Simple backrest, Brakes a little firmer, Post style hand grips for easier transfers.



## Client Portions as of July 2022:

Category 3 Rollator Base:	\$82.00
Category 2 Rollator Base:	\$63.50

Unilateral Hand Brake (Left or Right):	\$25.00
Forearm Trough Per Side:	\$50.00
Extra Wide (Greater than 18"):	\$25.00
Heavy Duty Package 1 (275-350lbs):	\$25.00
Heavy Duty Package 2 (350-400lbs):	\$37.50
Heavy Duty Package 3 (Greater than 400lbs):	\$50.00

\*Extra Width and Heavy Duty Package not combinable. ADP accepts whichever is greater.

### Accessories:

#### Nexus:

Nexus 3 Slow Down Brakes:	\$99.00
Nexus 3 Curb Climber:	\$12.95
Nexus Cup Holder:	\$25.95

#### Xpresso:

Xpresso Slow Down Brakes:	\$80.00
Xpresso Curb Climber:	\$19.95
Xpresso Wide/Low Slow Down Brakes:	\$96.00
Xpresso Cup Holder:	\$39.95

#### Trillium:

Trillium Cup Holder:	\$39.95
Trillium Slow Down Brakes:	\$79.95
Trillium Curb Climber:	\$19.95
Trillium 6 Slow Down Brakes:	\$79.95
Trillium 6 Curb Climber:	\$18.99

#### Nitro:

Nitro Cup Holder:	\$19.95
-------------------	---------

#### Stella:

Stella Slow Down Brakes:	\$82.00
Stella Footrests:	\$94.00
Stella Armrests:	\$61.00